



DOCTOR BROCHURE



+4822 675 24 94

grandortho@grandortho.pl



+420 730 136 494

patrik.goffa@grandortho.sk



+421 904 395 426

patrik.goffa@grandortho.sk



+37 066 040 050

monika.jakovleva@grandortho.lt



+38 067 352 93 98

grandortho@grandortho.com.ua

+38 067 362 93 98

Grand Ortho Sp. z o.o.
Bruszevska 28 str, 03-046 Warszawa



ORTHERO

CLEAREST WAY TO SMILE

ORTHERO



WORLDWIDE BRAND



10

YEARS



120+

EMPLOYEES



20+

COUNTRIES



4M+

ALIGNERS



60K+

HAPPY SMILES



ORTHERO

SUBMIT

Scan your patient and upload records to Orthero web portal

PLAN

Modify your treatment plan with sophisticated 3D Software

ORTHERO DIGITAL WORKFLOW

RETAIN

Long-lasting retention with HEROTAIN retainers


TREAT

Treat with high-tech material and achieve predictable results

MONITOR



Monitor your patient cooperation

SUBMIT




ORTHERO

PLAN



3D LAYER TECHNOLOGY



MADE
IN
GERMANY

PREDICTABLE



COMFORTABLE



STRONG



FLEXIBLE




PAINLESS



ORTHERO

PATIENT MONITORING APP



ORTHERO

HEROTAIN-ANTI RELAPSE SOLUTIONS

RELAPSE RECOVERY

POWERFUL

COMFORTABLE

POWERFUL

Herotain is thinner and %40 stronger, and does not wear out easily in long-term use

COMFORTABLE

Despite the strong structure, high elasticity makes it more comfortable to wear

RECOVERY

In case of a possible relapse, it has ability to recover up to 0.3mm

ORTHERO

ORTHERO

ORTHERO

DIFFERENCE



ULTRA-HIGH RESOLUTION PRINTING

MORE CLEAR



LESS BACTERIA



LESS COLORING



ORTHERO

DOUBLE LAYER ATTACHMENT TEMPLATE



EASY TO APPLY

Attachment templates inner thin foil is designed to minimize composite adhesion



SECOND CHANCE


Double-layer system allows a second chance to apply composite attachment

ORTHERO


OPTIMIZED CUTTING MARGIN

OPTIMIZED CUTTING MARGIN

Orthero takes advantage of the geometric stability on lingual-palatal region and having straight cut with gingival margin to have more active force and pressure. On the buccal-vestibular region, uses a scalloped cut without gingival margin to eliminate the geometric deviance occurrence after some stages of treatment.



CBCT INTEGRATION



Tooth Axis

Elimination of all assumption with real teeth axis set by CBCT integrated technology

Hinge Axis

Mandible rotation simulated by real hinge axis set from CBCT

Center of Resistance

Possibility to position attachments perfectly according to real center of resistance in terms of biomechanics

ORTHERO

AUTO-ROTATION



OPTIONS TO SET HINGE AXIS

MODJAW

SAM AXIOPRISA

CBCT WITH CONDYLES

CEPHALOMETRIC X-RAY




Our technology allows the auto-rotation of mandible from the real hinge axis set by various options.



Auto-Rotation from hinge axis obtains the most accurate result and prevents 1,5 to 2,5 mm sagittal axis deviations of virtual jump.

ORTHERO

PHASE BY PHASE



FLEXIBLE MANEUVER



If you need to change your treatment plan in the middle of the treatment, you can freely go another direction



SPLIT YOUR TREATMENT

You can divide your treatment in stages, and proceed phase by phase

ORTHERO

FAST STARTER KIT



ATTACHMENT TEMPLATE &
FIRST 2 STAGE OF ALIGNERS

**WORLD'S FASTEST
ALIGNER BRAND!**

ORTHERO



STANDART PACKAGES



AESTHETIC

Up to 8 Stages
1 Refinement
1 Year Warranty
1 Set HEROTAIN



FULL

Up to 16 Stages
2 Refinement
2 Years Warranty
1 Set HEROTAIN



ADVANCED LIGHT

Up to 25 Stages
3 Refinement
3 Year Warranty
1 Set HEROTAIN



ADVANCED

Unlimited Stage
Unlimited Refinement
5 Years Warranty
1 Set HEROTAIN

SPECIAL PACKAGES

EXPERT

Unlimited Stage
Unlimited Refinement
5 Year Warranty
Phase by Phase Option
CBCT Integrated Roots
Modjaw Integrated Hinge Axis



MINI

Unlimited Stage
Unlimited Refinement
2 Year Warranty
Phase by Phase Option
For age 6-12, in early mixed dentition
Permanent centrals, laterals and first molars must be erupted



ORTHERO

